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User Manual

BPB-M01H

Enjoy it.

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SPECIFICATIONS

- 1 Powerful Brushless High-Torque Motor With Quiet Glide Technology
- 2 Removable Battery Release
- 3 Rechargeable Lithium on Battery (2+ Hrs of Use Per Charge)
- 4 Multi-Speed Settings for Optimal Pressure
- 5 Cordless and Lightweight (2.5 Lbs) "Includes AC wall charger
- 6 ON/OFF Power Switch
- 7 Interchangeable Applicators(x4)

Fork

Bullet

Flat

Round

SPECIFICATIONS

Level1-20Hz	Level2-26.7Hz	Level3-33, 3H
Level4-40Hz	Level5-46, 7Hz	Level6-53, 3Hz

Charge:

100-240~50/60Hz 1.0A MAX

CHARGING

- Fully charge battery for up to six hours before first use either attached or detached from the device
- To charge, connect the dc end of the supplied 24V adapter to the battery's charging port and plug ac side cable into wall outlet
- LED band on battery lights pulsate to show battery level and indicate active charging
- With battery attached, LED lights on unit flash to show battery level and indicate active charging
- Full charge is indicated when all LED lights stop flashing and remain fully illuminated
- Battery may be recharged at any time and at any battery level
- It is not recommended to fully deplete the battery to the red led level
- The average run time is 2+ hours depending on speed level and pressure applied during use

OPERATING

- Align and attach battery into base of unit-ensure it is securely locked in place
- Power device on by moving the power switch on the bottom of the battery to the on position power device off by moving the power switch to the OFF position
- With power switch in ON position, press the power button once to turn device on to level one twice to turn device to level two and three times for level three. Pressing the power button again will turn the device off, the blue LED lights will indicate the corresponding power level
- With unit off, fully insert applicator by pressing in. Remove by pulling it on
- Massage desired body part at the desired speed applying pressure as needed
- To remove the battery, lightly press the battery release button with thumb firmly pull downwards on the battery handle with the other hand

MAINTENANCE CLEANING STORAGE

- Use a slightly damp towel to clean device and a dry with a soft doth
- For storage or travel, remove battery and store in supplied carton box

DISCLAIMER

Warnings & Caution

- Do not immerse in water. Keep away from liquids or heat sources. Keep ventilation ports free from dust and debris
- Do not remove screws or attempt to disassemble
- Do not operate continuously for more than one hour, allow device to rest for 30 minutes before re-using
- Unplug the charger from device battery after charging or prior to use
- For adult use only do not use if injured. Consult your doctor before using this product

SAFETY INSTRUCTIONS

DANGER: TO REDUCE RISKS OF ELECTRIC SHOCK, FIRE AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS

- For adult use only
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds per region
- Use the Massage Gun only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body
- Use only the applicators that provide the best desired effect
- Bruising can occur regardless of control setting or pressure applied check treated areas frequently and immediately stop at first sign of pain or discomfort
- Keep fingers hair or any other body part away from the shaft and backside of the applicator as pinching may occur
- Do not place any objects into the ventilation holes of the Massage Gun
- Do not immerse in water or allow water to enter the ventilation holes of the Massage Gun
- Do not drop or misuse the Massage Gun
- Only recharge with supplied 26 VDC Massage Gun charger
- Carefully examine device and battery before each use
- Do not tamper or alter Massage Gun in any way
- Never leave the Massage Gun operating or charging unattended

DISCLAIMER

PLEASE DO NOT USE THE MASSAGE GUN OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR FOR ANY OF THE FOLLOWING APPLIES

Pregnancy diabetes with complications such as neuropathy or retinal damage, wear of pace-makers, recent surgery epilepsy or migraines, herniated disks, spondylolisthesis, or spondylitis, recent joint replacement Or IUDs, metal pins or plates or any concerns about your physical health frail individuals and children should be accompanied by an adult when using any vibration device.

These contra indications do not mean that you are not able to use a vibration or exercise device but we advise you to consult a doctor first.

Ongoing research is being done to the effects of vibration exercise for specific medical disorders .This will likely result in the shortening of the list of contra indications as shown above, Practical experience has shown that integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advice of and in accompaniment with a doctor specialist or physiotherapist.

WARRANTY

IF WITHIN ONE YEAR FROM THE DATE OF PURCHASE, THIS PRODUCT FAILS DUE TO A DEFECT IN MATERIAL OR WORKMANSHIP WILL REPAIR OR REPLACE THE PRODUCT OR NECESSARY COMPONENTS, FREE OF CHARGE.

THIS WARRANTY EXCLUDES:

- (a) Damage caused by accident, abuse, mishandling or transport.
- (b) Sequence ton to unknown repair
- (c) Units not used in accordance with instructions
- (d) Damage exceeding the cost of the product.
- (e) Deterioration of the delivered product resulting from abnormal storage and or safeguarding conditions on the client's premises
- (f) Failure to provide the dated proof of purchase.

Some states and countries do not allow a limitation of damages. So the remaining limitation may not apply to you. This warranty guarantees specific legal rights, and but other rights may vary from country to country and from state to state. This warranty service is available through the operating policies and international distributors .

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