 **BLAUPUNKT**

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User Manual

BPB-M07HU

Enjoy it.

V1.0
007

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Contents

1. Important Safety Instruction	1
2. Product Feature	2
3. Product Parts	2
4. Massage Head	3
5. Installation	4
6. Operation Instruction & Charging	4-5
7. Specification	5
8. Cleaning and Maintenance	6
9. Trouble shooting	6
10.Storage	6

Important Safety Instruction

Please carefully read through the safety percussion before using the product.

Do not touch this part when starting the machine with massage head installed.

Do not touch this part when starting the machine without massage head installed.

- 1.Children need to be under supervision of adults to ensure that they do not play with the device.
- 2.This product is not intended for the users(including children) with undeveloped physical growth, sensing and mental capacity and those who is lack of experience or knowledge, unless it is safety used under the supervision or guidance of a person responsible for its safety.
- 3.Keep away from liquids or heat sources and keep the vents dust free.
- 4.This massage gun cannot be used during charging.
- 5.Please use this product according to the instructions.
- 6.Before use, please carefully check the product function according to the operation manual.

Product Features

1. Small, but strong power; Lightweight and compact; Easy to use by single hand
2. Touch screen button control, lightly touch to freely switch speed modes
3. 5 adjustable high percussion frequency up to 3500times/min
4. With patented silence and noiseless technology, just indulge yourself into quite and cozy relaxing hour
5. Standard brushless motor with powerful percussion to the deep fascia, to help muscle relaxation and deeply relieve muscle soreness
6. 3 massage heads replaceable, soothing muscles in various parts

Product Parts

Packing list: Mini Fascia Percussion Massage Gun

Accessories: massage headsx3, USB charge cablex1, satin bagx1
manual instruction(including warranty x1)

Massage Head

3 changeable massage heads equipped

Mini general massage head

Mini general massage head for better fascia stimulation suitable for multiple muscle groups!

General mode suitable for all muscle fascial parts.

Mini trigger points massage head

Mini trigger points massage head for deeper penetration into the “Trigger Points” and resolving these “adhesive pain points” one by one!

Dedicated Trigger points release mode for whole body trigger points relieving.

Mini achilles massage head

Mini achilles massage head fit adequately for ankles and achilles tendon, widely affected by exercisers and practitioners. Well performed in back muscles massage like hind neck and tendon.

Achilles tendon and ankle joint relaxation mode for relaxation before and after intense jumping and running exercise.

Installation

1. While installing, slightly rotate the massage head to load. (Figure1)
2. Turn the power switch to the “ON” position to turn on the device. (Figure2)
3. Press the screen touch button to switch between modes. (Figure3)
4. While pulling, slightly rotate the massage head to pull. Do not pull it out directly. (Figure4)

(Figure 1)

(Figure 2)

(Figure 3)

(Figure 4)

Operation Instruction & Charging

Instruction

Turn the power switch to the “ON” position, the device stands by
Speed mode: Short press the button to turn on, 1st speed by default; short press the button again to enter 2nd speed; short press the button again to enter 3rd speed; short press the button again to enter 4th speed; short press the button again to enter 5th speed.

Operation Instruction & Charging

Battery display

After power on, the battery power is < 10%, the red light is on;
the battery is 10%-25%, the light bar shows 2 grids, green light is always on;
the battery is 25%-50%, the light bar shows 3 grids, green light is always on;
the battery is 50%-75%, the light bar shows 4 grids, green light is always on;
the battery is 75%-100%, the light bar shows 5 grids, green light is always on.

Charging

1. Charge the battery before initial use.
(If the product is sufficiently charged, there is no need for charging).
2. When charging, insert the USB cable into the charger plug and connect it to the general household power outlet and connect the other end to the device charging port.
3. Battery level displays when charging:

Battery power 100%, the light bar shows 5 grids, green light is always on;

Battery power > 75%, the light bar shows 5 grids, the 1st, 2nd, 3rd, 4th grid are always on, the 5th grid is flashing;

Battery power 50%-75%, the light bar shows 4 grids, the 1st, 2nd, 3rd grid are always on, the 4th grid is flashing;

Battery power 25%-50%, the light bar shows 3 grids, the 1st and 2nd grid are always on, the 3rd grid is flashing;

Battery power 10%-25%, the light bar shows 2 grids, the 1st grid is always on, the 2nd grid is flashing;

Battery power less than 10%, the light bar shows 1 grid, red light flashing.

4. The device cannot be started up while charging.

Specification

Product Weight: about 0.59kg

5 Speeds adjustable

Rated Voltage: 5V-2A

Battery Volume: 2400mAh

Rated current: about 2000mA

Massage heads: 3 different massage heads

Battery Weight: around 143g

Charging time: around 4~6 hours

Discharging time: around 4~11 hours

Clean and Maintenance

Make Sure to switch off power and pull out the power plug before cleaning.

1. Clean up the surface of the device body. Wipe it with slightly wet cloth and dry with a soft cloth. Do Not clean with gasoline, cleaning naphtha, corrosive cleansing agents or solvents. Do Not flush (leach) with water.
2. Stored in a dry, well ventilated place with non-corrosive gas and no less than 250mm from the ground. The number of stacking layers does not exceed 5 layers.
3. If not in use for a long time, store in the package after cleaning and charge every six months.

Trouble shooting

No action

–Please check if the power switch is turned on.

Machine noise

–If there is foreign object or malfunction, please inspect it.

No screen display

–Please inspect the machine if there is malfunction

No charging

–Please check if the adaptor is plugged in properly and charging port is dirty.

Storage

1. When you go out or do not use the device for a long time, please switch off power to maintain product lifespan.
2. Please clean the device before storing to prevent mildew from non-using for a long time. Put on a plastic bag to prevent dust from entering.
3. Stored in the package provided when storing or traveling.
4. This product is suitable for car, train, ship and aircraft transport.
5. This product must be fixed during the transit to prevent severe vibration, collision, roll and fall and place it upward, lightly and somewhere rainproof.

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